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This manual is valid for the the ComboRehab-Series,
StimRehab-Series and UltraRehab-Series

This user manual is published by Shenzhen Dongdixin Technology Co., Ltd. and reserves the right to improve and amend it at any time without prior notice. Amendments may however be published in new editions of this manual.

All Rights Reserved.Rev.V1.1 © 2015, 20150529

Declaration of conformity

Shenzhen Dongdixin Technology Co.,LTD.declares that the ComboRehab-Series, StimRehab-Series and UltraRehab-Series complies with following normative documents:

**IEC60601-1,IEC60601-1-2,IEC60601-2-10,IEC60601-2-5,ISO7010
IEC61689,ISO14971,ISO10993-1,ISO10993-5,ISO10993-10**

**Complies with MDD 93/42/EEC and Amended by
directive 2007/47/EC requirements**

1. Foreword

1.1 Intended User/Operator

This manual has been written for the users of the Rehab-Series include ComboRehab-Series, StimRehab- Series and UltraRehab-Series. It contains general information on the operation, precautionary practices, and maintenance information. In order to maximize its use, efficiency, and the life of the system, please read this manual thoroughly and become familiar with the controls, as well as the accessories before operating the system. This device is designed to only be used by or under the supervision of persons using the medical device in the course of their work and in the framework of a professional health-care activity, who understand the benefits and limitations of electrotherapy and ultrasound therapy.

WARNING (USA only):

U.S.A. Federal Law restricts these devices to sale by, or on the order of, a physician or licensed practitioner. This device should be used only under the continued supervision of a physician or licensed practitioner.

Specifications put forth in this manual were in effect at the time of publication. However, owing to manufacturer her policy of continual improvement, changes to these specifications may be made at any time without obligation on the part of manufacturer.

2. Product Description

The Rehab-Series is a family of products for physical therapy, the Rehab-Series offers the practitioner a wide range of treatment options. The devices share an identical control panel equipped with a full colour touch panel means treatment setup has never been easier. A few simple key presses are all you need to quickstart a treatment. The User Interface intuitively groups and displays all the options for a modality setup on the large touch screen to ensure that treatment parameters can easily be selected and adjusted. The devices are mains powered and can optionally be equipped with a battery for mains independent operation. The family comprises the products described below.

StimRehab

The StimRehab is equipped with two or four completely identical electrotherapy channels. The electrotherapy channels can be used in combination (linked) or totally independent. A comprehensive set of current waveforms is available, targeting both pain management and muscle stimulation applications.

Protocols can run on linked or independent channels. With independent channels two or four different protocols can be performed simultaneously.

UltraRehab

The UltraRehab is an ultrasound therapy device. The device provides two positions for attachment of an ultrasound applicator. Depending on the device configuration ordered, the UltraRehab comes with an applicator with a large contact area, an applicator with a small contact area or with both applicators. The applicators can operate in continuous or pulsed mode at an ultrasound frequency of 1 MHz or 3 MHz. The lower frequency (1 MHz) penetrates deeper than a higher frequency (3 MHz), thus the practitioner can decide which frequency to use according to the condition and depth to be treated. Contact control suspends the application of ultrasonic energy when acoustical contact with the treatment area becomes insufficient. The applicators are suitable for subaqueal treatments.

ComboRehab

The ComboRehab is a combination device, combining the functions of the StimRehab and the UltraRehab in a single device. With the ComboRehab the simultaneous application of ultrasound and electrotherapy (combination therapy) is also possible. The remaining electrotherapy channel can then be used independently.

VAM200

Electrotherapy can be applied through standard or vacuum electrodes. With vacuum electrodes the VAM200 generates the vacuum through which the vacuum electrodes are attached to the patient. The device is placed beneath the StimRehab or ComboRehab, from which its power is derived and through which it is also operated.

MTM200

The MTM200 is a two channel electrotherapy module intended to upgrade the StimRe

hab and ComboRehab to Four Channel Electrotherapy or Combination Therapy Systems. This module is designed for use with the StimRehab or ComboRehab only.

EMG200

The EMG200 is designed for use with the StimRehab and the ComboRehab. With EMG electrodes the EMG200 generates the sEMG (Surface Electromyography) , ETS (Surface Electromyography with Triggered Stimulation) through which the EMG electrodes are attached to the patient.

BTM200

The BTM200 Battery Module is designed for use with the Rehab-Series systems to create a battery powered Therapy System. No additional Software is required for the Module as the System automatically recognizes its presence and activates all necessary software inherent in the System.

Table 1 Difference for the models in Rehab-Series

Name	Model	Mainframe	Vacuum module (VAM200)	EMG module (EMG200)	Battery module (BTM200)
ComboRehab ²	CT2200 (two channels electrotherapy, ultrasound and combo)	Mainframe CT2200			
ComboRehab ² Vac	CT2201 (two channels electrotherapy, ultrasound and combo+vacuum module)	Mainframe CT2200	✓		
ComboRehab ² Bio	CT2202 (two channels electrotherapy, ultrasound and combo+EMG module)	Mainframe CT2200		✓	
ComboRehab ² Plus	CT2203 (two channels electrotherapy, ultrasound and combo+battery module)	Mainframe CT2200			✓

Name	Model	Mainframe	Vacuum module (VAM200)	EMG module (EMG200)	Battery module (BTM200)
ComboRehab² VB	CT2204 (two channels electrotherapy, ultrasound and combo+vacuum module +EMG module)	Mainframe CT2200	✓	✓	
ComboRehab² Vac Plus	CT2205 (two channels electrotherapy, ultrasound and combo+vacuum module+ battery module)	Mainframe CT2200	✓		✓
ComboRehab² Bio Plus	CT2206 (two channels electrotherapy, ultrasound and combo+battery module + EMG module)	Mainframe CT2200		✓	✓
ComboRehab⁴	CT2400 (four channels electrotherapy, ultrasound and combo)	Mainframe CT2400 (CT2200 +MTM200)			
ComboRehab⁴ Vac	CT2401 (four channels electrotherapy, ultrasound and combo+vacuum module)	Mainframe CT2400 (CT2200 +MTM200)	✓		
ComboRehab⁴ Bio	CT2402 (four channels electrotherapy, ultrasound and combo+EMG module)	Mainframe CT2400 (CT2200 +MTM200)		✓	
ComboRehab⁴ Plus	CT2403 (four channels electrotherapy, ultrasound and combo+battery module)	Mainframe CT2400 (CT2200 +MTM200)			✓
StimRehab²	MT2000 (two channels electrotherapy)	Mainframe MT2200			

Name	Model	Mainframe	Vacuum module (VAM200)	EMG module (EMG200)	Battery module (BTM200)
Stim Rehab ² Vac	MT2201 (two channels electrotherapy + vacuum module)	Mainframe MT2200	✓		
Stim Rehab ² Bio	MT2202 (two channels electrotherapy + EMG module)	Mainframe MT2200		✓	
Stim Rehab ² Plus	MT2203 (two channels electrotherapy + battery module)	Mainframe MT2200			✓
Stim Rehab ² VB	MT2204 (two channels electrotherapy + vacuum module + EMG module)	Mainframe MT2200	✓	✓	
Stim Rehab ² Vac	MT2205 (two channels electrotherapy + vacuum module + battery module)	Mainframe MT2200	✓		✓
Stim Rehab ² Bio Plus	MT2206 (two channels electrotherapy + battery module + EMG module)	Mainframe MT2200		✓	✓
Stim Rehab ⁴	MT2400 (four channels electrotherapy)	Mainframe MT2400 (MT2200 + MTM200)			
Stim Rehab ⁴ Vac	MT2401 (four channels electrotherapy + vacuum module)	Mainframe MT2400 (MT2200 + MTM200)	✓		

Name	Model	Mainframe	Vacuum module (VAM200)	EMG module (EMG200)	Battery module (BTM200)
Stim Rehab ⁴ Bio	MT2402 (four channels electrotherapy +EMG module)	Mainframe MT2400 (MT2200 +MTM200)		✓	
Stim Rehab ⁴ Plus	MT2403 (four channels electrotherapy +battery module)	Mainframe MT2400 (MT2200 +MTM200)			✓
UltraRehab	UT2200 (two channels ultrasound)	Mainframe UT2200			
UltraRehab Plus	UT2201 (two channels ultrasound +battery module)	Mainframe UT2200			✓

Note:Remarks: "✓" means that the device is equipped with these module.

3. Precautionary Instructions

In this section general Warnings and Precautions are listed, that you should be aware of when using the Rehab-Series. See also chapter 4.1 for Warnings and Precautions that are application specific.

WARNING:

- Federal law (USA only) restricts this device to sale by, or on the order of, a physician or licensed practitioner. This device should be used only under the continued supervision of a physician or licensed practitioner.
- Make certain that the unit is electrically grounded by connecting only to a grounded electrical service receptacle conforming to the applicable national and local electrical codes.
- Do not operate the unit in an environment of short-wave or micro-wave diathermy.
- The device is designed for indoor use only. It is prohibited to use the device in a location where explosion or water intrusion risk are present and in dusty or humid environment. It is prohibited to use the device in spaces where flammable anaesthetics oxidizing gases (O₂, N₂O) and other flammable gases or vapors are present.
- This device should be kept out of the reach of children.
- Place the device out of direct sunlight and strong electromagnetic fields of surrounding devices (diathermy, X-rays, mobile phones and other radio-frequency equipment) to prevent unwanted interference. If unwanted interference occurs, place the device farther from the source of interference or contact the Nu-Tek authorized service.
- Before administering any treatment to a patient you should become acquainted with the operating procedures for each mode of treatment available, as well as the indications, contraindications, warnings and precautions. Consult other resources for additional information regarding the application of electrotherapy and Ultrasound.
-  No modification of this equipment is allowed! Do not try to open or remove the protective covers or disassembly the device for any reason. There is a danger of electric shock and serious injury. All service actions must be done by an authorized Nu-Tek service only; otherwise Nu-Tek bear no responsibility for further operation of the device.

- Never use the accessories connector and other connectors to plug in anything else except transducers and cables sold by the manufacturer as replacement parts for internal components. There is a serious risk of electric shock and serious damage to the device!
- This device is not designed to be use in an MRI Environment and should be removed prior to MRI exposure.

 **CAUTION:**

- Keep yourself informed of the contraindications.
- Read the User's Manual carefully and become familiar with all its safety requirements, operating procedures and maintenance instruction prior to use of the device. It is prohibited to use the device and its accessory in any manner that is not in accordance with the User's Manual. Know the limitations and hazards associated with using any electrical stimulation device. Observe the precautionary and operational stickers placed on the unit.
- DO NOT use sharp objects such as a pencil point or ballpoint pen to operate the buttons on the control panel.
- Before each therapy check carefully the device and its accessories (cables, connectors, electrodes, ultrasound heads, controls, touch screen) for any mechanical, functional or other damage. If any faults or anomalies in the device function are found, stop using the device immediately .
- Handle the ultrasound applicator with care. Inappropriate handling of the ultrasound applicator may adversely affect its characteristics.
- Inspect the ultrasound applicator for cracks which may allow the ingress of conductive fluid before each use.
- This unit should be operated in temperatures between 10 °C and 40 °C (50 °F and 104 °F), with a Relative Humidity ranging from 30%-85% non condensing.
- Do not expose the unit to direct sunlight, heat radiated from a heat radiator, excessive amounts of dust, moisture, vibrations and mechanical shocks.
-  The device has applied parts of the BF (Body Floating) type-i.e. parts which come into direct physical contact with the patient during normal device use. This includes the electrodes for electrotherapy and applicators for ultrasound therapy.
- The device heats up during operation and therefore must not be located near devices that heat up or produce heat. The device is cooled by forced air circulation. The cooling vents are located on the rear and side panel of the device and they must not be covered. When placing the device, leave at least 10 cm of space behind the rear panel.

- It is prohibited to place any objects that produce heat or objects that contain water or other liquid on the device. If in the case of ingress of liquids, unplug the unit from the mains supply and have it checked by an authorized person (see the paragraph on maintenance).
- Before administering any treatment to a patient you should become acquainted with the operating procedures for each mode of treatment available, as well as the indications, contra-indications, warnings and precautions. Consult other resources for additional information regarding the application of electrotherapy and ultrasound therapy.
- The device is designed to only be used by or under the supervision of persons using the medical device in the course of their work and in the framework of a professional healthcare activity, who understand the benefits and limitations of electrotherapy and ultrasound therapy, i.e. "professional users".

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9.1 Cleaning

9.1.1 Cleaning of the device

Switch off the device and disconnect it from the power supply. The apparatus can be cleaned with a damp cloth. Use lukewarm water and a non-abrasive liquid household cleaner (no abrasive, no alcohol content solution). If a more sterile cleaning is needed, use a cloth moistened with an antimicrobial cleaner.

CAUTION:

Do not submerge the apparatus in liquids. Should the unit accidentally become submerged, contact the dealer or Authorized Service center immediately. Do not attempt to use a system that has been wet inside until inspected and tested by a Service Technician Certified by Authorized Service center. Do not allow liquids to enter the ventilation holes.

9.1.2 Cleaning of display panel

Use a soft and dry cotton cloth or micro fiber tissue to clean the panel. To remove fingerprints or grease, use a non-abrasive glass cleaning agent. Apply a small amount of the cleaning agent to a soft cotton cloth and then carefully clean the panel.

CAUTION:

- Do not spray the cleaning agent directly on the glass panel.
- Do not use cleaning agents that contain strong alkalis, lye, acid, detergents with fluoride or detergents with ammonia.

9.1.3 Cleaning the electrodes

- Apply the protective backing to the tacky side of the electrode. Place the electrode on the side of the protective backing that is labeled with the word on.
- It may be helpful to improve repeated application by spreading a few drops of cold water over the adhesive and turn the surface up to air dry. Over Saturation with water will reduce the adhesive properties.
- Between uses, store the electrodes in the reusable bag in a cool dry place.

CAUTION:

- The electrodes are intended for single patient use only .
- If irritation occurs, discontinue use and consult your clinician.
- Always use the electrodes with CE mark, or are legally marketed in the US under 510(K) procedure.

9.1.4 Cleaning the lead wires and cables

Periodically wipe the lead wires clean with a cloth dampened in a mild soap solution, and then gently wipe them dry . Use of rubbing alcohol on the lead wires will damage the insulation and dramatically shorten their life.

9.1.5 Ultrasound applicator

To prevent corrosion, clean and dry the contact surface immediately after use. Make sure that no ultrasound gel remains on the applicator. We further recommend cleaning the applicator and cable daily, using lukewarm water. The applicator can be disinfected using a cloth moistened with a 70% alcohol solution. Check the applicator and cable regularly for damage.

9.1.6 Vacuum electrodes and sponges

The vacuum electrodes and sponges should be cleaned with lukewarm water. In the case of persistent dirt, and for disinfection, a 70% alcohol solution may be used.

Sponges should be replaced regularly. It is recommended to keep sponges and a spare electrode in stock.

Calcium scale can be deposited on the metal surfaces of the electrodes. This has an insulating effect. In order to maintain optimum conductivity, these surfaces should be regularly cleaned and polished.

9.1.7 Vacuum cables

Clean the vacuum cable with a damp cloth. Use lukewarm water and a non-abrasive household cleaning agent. Do not use an alcohol solution. Check the cable regularly for damages and/or bad electrical contact. We advise, keeping a spare vacuum cable in stock.

9.1.8 Cleaning Vaginal/Rectal probes

Carefully clean the probe after use. Wash the probe gently in mild soapy water, rinse and make sure the probe is completely dry before returning to storage in the plastic bag.

9.1.9 Cleaning the water reservoir and hoses

- Detach the vacuum cups from the vacuum cables.
- Place a container filled with a cleaning liquid below the system.
- Place the peripheral ends of the cables in the container.
- Go to System Settings and select Tank Cleaning.
- The water reservoir will be filled with the cleaning liquid until the water reservoir is full.

9.2 Warning Messages, Error Messages

9.2.1 Prompt tone

error code	reason	Content		System processing
		Tooltip	Voice	
001	Ultrasound applicar is without load	N/A	Prompt tone	Output pause, system recover output when detecting load.
002	Parameters is set upper limit	N/A	Prompt tone	Parameters remain in upper limit
003	Parameters is set lower limit	N/A	Prompt tone	Parameters remain in lower limit
004	Choose occupied channel operation	N/A	Prompt tone	The operation is invalid
005	Touch operation that is not allowed	N/A	Prompt tone	The operation is invalid
006	Touch option that cannot be changed	N/A	Prompt tone	The operation is invalid

9.2.2 Warning tone

error code	reason	Content		System processing
		Tooltip	Voice	
101	Electrical stimulation channel 1 is without load	Bad contact quality on channel 1.Check pads and lead wires.	Warning tone	System Stop Electrical stimulation output
102	Electrical stimulation channel 2 is without load	Bad contact quality on channel 2.Check pads and lead wires.	Warning tone	System Stop Electrical stimulation output
103	Electrical stimulation channel 3 is without load	Bad contact quality on channel 3.Check pads and lead wires.	Warning tone	System Stop Electrical stimulation output
104	Electrical stimulation channel 4 is without load	Bad contact quality on channel 4.Check pads and lead wires.	Warning tone	System Stop Electrical stimulation output
105	The program saved is more than system capacity	Attempting to save Favorite Protocols after system memory has reached the maximum allowed	Warning tone	The operation is invalid
106	System detect Ultrasound handle is not connect in Ultrasound treatment	Ultrasound Applicator disconnected from system during treatment session	Warning tone	System Stop Ultrasound treatment output

107	System detect Ultrasound handle is not connect when the Ultrasound treatment is selected.	Attempting to perform Ultrasound treatment with no Applicator connected to the system	Warning tone	System Stop Ultrasound treatment output
108	Ultrasound applicar overheat	Ultrasound Applicator is too hot.	Warning tone	System Stop Ultrasound treatment output
109	Electrical stimulation channel 1 short circuit	Overcurrent on channel 1. Check pads and lead wires.	Warning tone	System Stop Electrical stimulation output
110	Electrical stimulation channel 2 short circuit	Overcurrent on channel 2. Check pads and lead wires.	Warning tone	System Stop Electrical stimulation output
111	Electrical stimulation channel 3 short circuit	Overcurrent on channel 3. Check pads and lead wires.	Warning tone	System Stop Electrical stimulation output
112	Electrical stimulation channel 4 short circuit	Overcurrent on channel 4. Check pads and lead wires.	Warning tone	System Stop Electrical stimulation output
113	Battery level is Low	Battery level too Low	Warning tone	The operation is invalid
114	Battery is not enough	The battery is insufficiently charged to complete the treatment at the currently set therapy levels.	Warning tone	The operation is invalid
115	water reservoir is full	The water separation tank of the Vacotron is full.	Warning tone	The operation is invalid

116	Vacuum adsorption module e leak	There probably is a leak in the vacuum system.	Warning tone	The operation is invalid
117	Device inside Overtemperature	System is too hot.	Warning tone	System stop output
118	Battery overtemperature	Battery is too hot.	Warning tone	System stop output

9.3 Maintenance

9.3.1 User Maintenance

9.3.1.1 Technical Maintenance

On request a service manual can be made available containing: spare part list, descriptions, calibration instructions and other information which will assist the user’s qualified technical personnel to repair those parts of the equipment which are designated by the manufacturer as repairable.

⚠ CAUTION:

- Electrical safety of the device relies on a properly earthed electrical connection via the power cord. It is therefore necessary to have this connection checked annually.
- To ensure continued compliance with the 21 CFR 1050.10 standard, this unit should be adjusted and safety tested once each year. Procedures laid down in the service manual should be followed. This may be carried out by your supplier, or by another agency, authorized by the manufacturer. It is also recommended that a service history record is maintained. In some countries this is even obligatory.
Use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous exposure to ultrasonic energy.

⚠ WARNING:

- This unit operates with high voltages. No attempt should be made to disassemble the unit. Maintenance and repair should be carried out by authorized personnel only. The manufacturer will not be held responsible for the results of maintenance or repairs by unauthorized persons.

All other technical maintenance is restricted to authorized Nu-Tek maintenance personnel. Authorized service personnel can make use of 1498770 Service manual REHAB SERIES

9.4 Troubleshooting

- Replace lead wires annually.
- Please follow the directions on the electrode packaging for the care of electrodes. The life of the electrodes varies, depending on skin conditions, skin preparation, storage and climate. Replace electrodes.
- that no longer stick.
- **NOTE:** If the following measures fail to alleviate the problem, please call the authorized agency or your supplier

Problem	Possible cause	Solution
Displays fail to light up	Adapter contact failure	Ensure adapter is connect. Check the following contacts: <ul style="list-style-type: none"> • All contacts are in place. • All contacts are not broken. • Ensure that adapter is Connected.
Displays fail to light up	Electrodes 1. Dried out or contaminated 2. Placement	1. Replace. 2. Electrodes must be a minimum of 2 inches apart.
	Lead wires Old/worn/damaged	Lead wires Old/worn/damaged
Stimulation stops	Poor electrode contact	Reapply electrodes, secure firmly
	Damaged or worn electrodes or lead wires	Replace

Stimulation is uncomfortable.	Intensity is too high	Decrease intensity .
	Electrodes are too close together	Reposition the electrodes.
		Electrodes must be a minimum of 2 inches apart.
	Damaged or worn electrodes or lead wires	Replace
	Electrode active area size is too small	Replace electrodes with ones that have an active area no less than 25.0cm .
Stimulation is ineffective	Improper electrode	Reposition electrode
	Unknown	Contact clinician

9.5 End of life

The REHAB SERIES contains materials that can be recycled and/or are noxious to the environment. Specialized companies can dismantle the unit and sort out these materials. When you dispose of the unit, find out about local regulations concerning waste management

10. Specifications

10.1 Ultrasound parameters

Frequency.....	1MHZ, ± 10%; 3 MHZ, ± 10%
Duty Cycles.....	10%-100%, Continues, stepping 10%
Pulse duration.....	1 – 56 ms ± 10 % (set by duty cycle)
PulseFrequency.....	16, 48, 100Hz

Output Power

Duty factor ≥ 80% for 5cm ²	0.5W-10.0W
Duty factor ≤ 70% for 5cm ²	0.5W-15.0W
Duty factor ≥ 80% for 1cm ²	0.1W-2.0W
Duty factor ≤ 70% for 1cm ²	0.1W-3.0W
Output accuracy.....	± 20% (for any level above 10% of maximum)

Peak out Amplitude Duty

Duty factor ≥ 80%.....	2.0W/cm ²
Duty factor ≤ 70%.....	3.0W/cm ²
Treatment timer.....	0 - 30 min ± 0.1 min

5 cm² Applicator

ERA (Effective Radiation Area)	5cm ²
Beam type:	
1 MHz.....	Collimating
3 MHz.....	Collimating
BNR (Beam Non-uniformity Ratio).....	5:1 maximum

1 cm² Applicator

ERA (Effective Radiation Area).....	1cm ²
Beam type:	
1 MHz.....	1cm ²
3 MHz.....	1cm ²
BNR (Beam Non-uniformity Ratio).....	5:1 maximum

10.2 Stimulator output parameters

10.2.1 IFC-4P: IFC(Interferential) Traditional (4 Pole)

Interferential Current is a medium frequency waveform. Current is distributed through of two channels (four electrodes). The currents cross each other in the body at the area requiring treatment. The two currents interfere with each other at this crossing point, resulting in a modulation of the intensity (the current intensity increases and decreases at a regular frequency).

Vector Scan.....Auto:20%-100%,stepping 20%;manual:0°-90°,stepping 15°
 Output Mode.....Electrodes
 Carrier Frequency.....2-10KHz, stepping 0.5KHz
 Beat High.....(Beat L.)-200 Hz,stepping 1Hz
 Beat Low.....1-(Beat H.)Hz,stepping 1Hz
 Intensity.....CC:0-100mA,stepping 0.5mA;CV:0-100V,stepping 0.5V
 Treatment Time.....1-60 minutes

10.2.2 IF-2P: IFC(Interferential) Premodulated (2 Pole)

Premodulated Current is a medium frequency waveform. Current comes out of one channel (two electrodes). The current intensity is modulated: it increases and decreases at a regular frequency (the Amplitude Modulation Frequency).

Carrier Frequency.....2-10KHz, stepping 0.5KHz
 Beat High.....(Beat L.)-200 Hz,stepping 1Hz
 Beat Low.....1-(Beat H.)Hz,stepping 1Hz
 Intensity.....CC:0-100mA,stepping 0.5mA;CV:0-100V,stepping 0.5V
 Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Treatment Time.....1-60 minutes
 Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1, 2/2,5/5,8/8,custom:0-9

10.2.3 Biphasic (TENS)

The Asymmetrical Biphasic and the Symmetrical Biphasic waveform are often used in TENS (Transcutaneous Electrical Nerve Stimulation) applications. The TENS has a short pulse duration. It is capable of strong stimulation of the nerve fibers in the skin as well as of muscle tissue. Because of its short pulse, the patient typically tolerates the current well, even at relatively high intensities. The Alternating Rectangular waveform is an interrupted biphasic current with a rectangular pulse shape. This waveform is commonly used as a pain management application.

10.2.3.1 TENS Asymmetrical

Frequency1-250Hz, stepping 1Hz
 F.M.0-250Hz, stepping 1Hz,cycle 16s
 P.Dur.20 μ s-1000 μ s,stepping 5 μ s
 Intensity.....CC:0-200mA,stepping 0.5mA;CV:0-100V,stepping 0.5V
 Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Treatment Time.....1-60 minutes
 Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1, 2/2,5/5,8/8,custom:0-9

10.2.3.2 TENS Asymmetrical Burst

Frequency1-250Hz, stepping 1Hz
 Burst.1-9bps, stepping 1bps
 P.Dur.20 μ s-1000 μ s,stepping 5 μ s
 Intensity.....CC:0-200mA,stepping 0.5mA;CV:0-100V,stepping 0.5V
 Treatment Time.....1-60 minutes
 Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1, 2/2,5/5,8/8,custom:0-9

10.2.3.3 TENS Symmetrical

Frequency1-250Hz, stepping 1Hz
 F.M.0-250Hz, stepping 1Hz
 P.Dur.20 μ s-1000 μ s,stepping 5 μ s
 Intensity.....CC:0-200mA,stepping 0.5mA;CV:0-100V,stepping 0.5V

Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Treatment Time..... 1-60 minutes
 Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1, 2/2,5/5,8/8,custom:0-9

10.2.3.4 TENS Symmetrical Burst

Frequency.....1-250Hz, stepping 1Hz
 Burst 1-9bps, stepping 1bps
 P.Dur.....20 μ s-1000 μ s,stepping 5 μ s
 Intensity.....CC:0-200mA,stepping 0.5mA;CV:0-100V,stepping 0.5V
 Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Treatment Time..... 1-60 minutes
 Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1, 2/2,5/5,8/8,custom:0-9

10.2.3.4 TENS Alternating Rec.

Frequency1-250Hz, stepping 1Hz
 F.M.....0-250Hz, stepping 1Hz
 P.Dur.....20 μ s-1000 μ s,stepping 5 μ s
 Intensity.....CC:0-200mA,stepping 0.5mA;CV:0-200V,stepping 0.5V
 Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Treatment Time..... 1-60 minutes
 Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1, 2/2,5/5,8/8,custom:0-9

10.2.3.5 TENS Alternating Rec.Burst

Frequency1-250Hz, stepping 1Hz
 Burst..... 1-9bps, stepping 1bps
 P.Dur.....20 μ s-1000 μ s,stepping 5 μ s
 Intensity.....CC:0-200mA,stepping 0.5mA;CV:0-200V,stepping 0.5V
 Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Treatment Time..... 1-60 minutes
 Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1, 2/2,5/5,8/8,custom:0-9

10.2.4 Russian

Russian Current is a rectangle waveform, delivered in bursts or series of pulses. This method was claimed by its author (Kots) to produce maximal muscle strengthening effects without significant discomfort to the patient.

Carrier Freq.....2-10KHz
 Burst Frequency.....20-100Hz,stepping 5Hz
 Duty cycle.....10%-50%,stepping 10%
 Intensity.....CC:0-200mA,stepping 0.5mA;CV:0-200V,stepping 0.5V
 Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Treatment Time.....1-60 minutes
 Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1, 2/2,5/5,8/8,custom:0-9

10.2.5 Microcurrent

Microcurrent is a monophasic waveform of very low intensity. The literature reports beneficial effects of this waveform in the treatment of wounds. The physiological working mechanism of this effect is as yet not clearly understood. It is thought to stimulate tissue healing by stimulating the 'current of injury', a current which naturally occurs in healing tissue.

Frequency.....0.1-1000Hz,stepping 0.1Hz/1Hz
 Polarity.....Positive,negative
 Intensity.....CC:0-1000 μ A,stepping 5 μ A
 Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Treatment Time.....1-60 minutes
 Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1,2/2,5/5,8/8,custom:0-9

10.2.6 Faradic

10.2.6.1 Trabert

It is a monophasic waveform with a phase duration of 2 ms and a pause of 5 ms resulting in a frequency of approximately 143 Hz.

Polarity.....Positive,negative

Phase Duration.....Positive,negative
 Phase Duration.....2ms
 Interval.....5ms
 Intensity.....CC:0-70mA,stepping 0.5mA;CV:0-70V,stepping 0.5V
 Treatment Time.....1-60 minutes

10.2.6.2 Rectanglar

Frequency.....0.2Hz-200Hz.stepping 0.1Hz/1Hz
 Polarity.....Positive,negative
 P.Dur.....100 μ s-1000ms,steping 100 μ s-100ms
 Intensity.....CC:0-80mA,stepping 0.5mA;CV:0-80V,stepping 0.5V
 Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Treatment Time.....1-60 minutes
 Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1, 2/2,5/5,8/8,custom:0-9

10.2.6.3 Triangular

Frequency.....0.2Hz-200Hz.stepping 0.1Hz/1Hz
 Polarity.....Positive,negative
 P.Dur.....100 μ s-1000ms,steping 100 μ s-100ms
 Intensity.....CC:0-80mA,stepping 0.5mA;CV:0-80V,stepping 0.5V
 Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Treatment Time.....1-60 minutes
 Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1, 2/2,5/5,8/8,custom:0-9

10.2.7 Diadynamic

The Diadynamic waveforms are rectified alternating currents. The alternating current is modified (rectified) to allow the current to flow in one direction only.

Mode.....MF,dF,LP,CP,CP.id
 Polarity.....Positive,negative
 Intensity.....CC:0-70mA,stepping 0.5mA; CV:0-70V,stepping 0.5V
 Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Treatment Time.....1-60 minutes

Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1, 2/2,5/5,8/8,custom:0-9

10.2.8 High Voltage

The High Voltage Pulsed Current (HVPC) has a very brief pulse duration characterized by two distinct peaks delivered at high voltage. The waveform is monophasic (current flows in one direction only). The high voltage causes a decreased skin resistance making the current comfortable and easy to tolerate.

Frequency.....1-120Hz,stepping 1Hz
 Polarity.....Positive,negative
 P.Dur.....100 μ s
 Intensity.....CV:0-500V,stepping 5V
 Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Treatment Time.....1-60 minutes
 Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1, 2/2,5/5,8/8,custom:0-9

10.2.9 NMS

10.2.9.1 NMS

Mode.....NMS,NMS Burst
 Frequency.....1-250Hz,stepping 1Hz
 P.Dur.....20 μ s-400 μ s,stepping 5 μ s
 Intensity.....CC:0-200mA,stepping 0.5mA;CV:0-200V,stepping 0.5V
 Channel mode.....Single channel Mode,Reciprocal Mode,Co-Contraction Mode
 Treatment Time.....1-60 minutes
 Cycle Time.....continues,10/10,10/20,10/30,10/50,custom:0-60
 Ramp.....0/0,1/1, 2/2,5/5,8/8,custom:0-9

10.2.10 Galvanic

10.2.10.1 Continuous

Polarity.....Positive,negative
 Intensity.....CC:0-80mA,stepping 0.5mA;CV:0-80V,stepping0.5V

10.2.10.2 Interrupted

Frequency.....8KHz fixed
 Polarity.....Positive,negative
 Duty cycle.....90% fixed
 Intensity.....CC:0-80mA,stepping 0.5mA;CV:0-80V,stepping0.5V

10.3 Parameter Limit

For security, the device has some some limited for electrotherapy.The maximum intensity has a relationship with the frequency and pulse duration following below table:

Pulse frequency	phase duration	Max Current output
<100	<300 μ S	200mA
<100	300-500 μ Sec	150mA
<100	500-1000 μ Sec	100mA
100-250	<500 μ Sec	100mA
100-250	500-1000 μ Sec	70mA

10.4 Technical Data

Power supply.....100V-240V, 47Hz-63Hz, 1.35A
 Power output.....15V , 4A Max
 Dimensions.....256x180x124mm(LxWxH)
 Operating Environmental:
 Temperature.....10°C(50°F) to 40°C(104°F)
 Relative humidity.....30%-85%
 Atmosphere pressure.....700hPa-1060hPa
 Storage Environmental:
 Temperature.....-10°C(14°F) to 55°C(131°F)
 Relative humidity.....10%-90%
 Atmosphere pressure.....700hPa-1060hPa

10.5 Safety and Performance standards

IEC 60601-1

General requirements for the safety of electrical medical systems, including Annex 1, national differences for Australia, Canada and the United States.

Safety class according to IEC 60601-1

class I type BF 

IEC 60601-2-5

Particular requirements for the safety of ultrasonic therapy equipment.

IEC 60601-2-10

Particular requirements for the safety of nerve and muscle stimulators.

This equipment complies with all requirements of the Medical Device Directive (93/42/EEC).

Medical device classification

IIa

21 CFR 1050.10

This equipment complies with all requirements of 21 CFR 1050.10, Performance Standard for Ultrasonic Therapy devices.

21 CFR 898

This equipment complies with all requirements of 21 CFR 898, Performance Standard for electrode lead wires and patient leads.

10.6 EMC details

Medical electrical devices such as the Rehab-series are subject to special precautions with regard to electromagnetic compatibility (EMC) and must be installed and commissioned in accordance with the EMC advice given in the instructions for use and accompanying documents.

Portable and mobile RF communication systems (e.g. mobile phones) may interfere with medical electrical Rehab-series.

The Rehab-series should only be operated with the original mains cable specified in the list of contents delivered.

Operating the device with any other mains cable can lead to increased emissions or reduced interference immunity of the device.

Guidelines and manufacturer's declaration – electromagnetic interference		
The REHAB SERIES device is intended for operation in an electromagnetic environment as indicated below. The customer or user of the REHAB SERIES unit should ensure that it is operated in such an environment.		
Interference tests	Conformity	Electromagnetic environment guideline
RF emissions according to CISPR 11	Group 1	The Rehab-series device uses RF energy solely for its internal functioning. Its RF emission is therefore very low and it is unlikely that this will cause interference to neighbouring electronic Rehab-series
RF emissions according to CISPR 11	Class A	The Rehab-series device is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions according to IEC 61000-3-2	Class A	
Voltage fluctuation emissions and flicker according to IEC 61000-3-3	Conforms	

The device should not be used when placed immediately next to or stacked on top of other devices. If operation is necessary when immediately next to or stacked on top of other devices, the device should be monitored to ensure it is operating as intended in this arrangement.

Guidance and manufacturer's declaration – Electromagnetic immunity			
The REHAB SERIES device is intended for use in the electromagnetic environment specified below. The customer or the user of the REHAB SERIES device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - Guidance
Electrostatic discharge (ESD) to IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst to IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input /output lines	± 2 kV for power supply lines not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 6100-4-5	± 1 kV differential mode ± 2 kV common mode	± 1 kV differential mode ± 2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT for 0.5 cycle) 40% UT (60% dip in UT for 5 cycles) 70% UT (30% dip in UT for 25 cycles) <5% UT (>95% dip in UT for 5 seconds)	<5% UT (>95% dip in UT for 0.5 cycle) 40% UT (60% dip in UT for 5 cycles) 70% UT (30% dip in UT for 25 cycles) <5% UT (>95% dip in UT for 5 seconds)	Mains power quality should be that of a typical commercial or hospital environment. If the use of the REHAB SERIES device requires continued operation during mains power interruptions, it is recommended to install a battery.
Power frequency (50/60 Hz) magnetic field to IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: UT is the AC mains voltage prior to application of the test level.			

The main features of the REHAB SERIES devices are as follows: interference-free delivery of shockwaves, interference-free control of all functions. Uninterrupted operation is not required with the use intended.

Guidelines and manufacturer's declaration – electromagnetic interference immunity			
The REHAB SERIES device is intended for operation in the electromagnetic environment specified below. The customer or user of the REHAB SERIES should ensure that it is used in such an environment.			
Interference immunity tests	IEC 60601-test level	Compliance level	Electromagnetic environment-guidelines
<p>Conducted RF disturbance variables according to IEC 61000-4-6</p> <p>Radiated RF disturbance variables according to IEC 61000-4-3</p>	<p>3 Veffective value 150 kHz to 80MHz</p> <p>3 V/m 80 MHz to 2.5GHz</p>	<p>3 Veffective value 150 kHz to 80 MHz</p> <p>3 V/m 80 MHz to 2.5 GHz</p>	<p>Portable and mobile radio should not be used any closer to the REHAB SERIES devices, including cables, than the recommended separation distance calculated from the equation applicable to the transmission frequency.</p> <p>Recommended separation distance: $d = 1.2 \sqrt{P}$ $d = 0.35 \sqrt{P}$ for 80 MHz to 800 MHz $d = 0,7 \sqrt{P}$ for 800 MHz to 2.5 GHz</p> <p>Where P is the rated power of the transmitter in Watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>According to an investigation in situa, the field strength of stationary radio transmitters should be less than the compliance level at all frequencies.</p> <p>Interference may occur in the vicin-</p> 

NOTE 1 At 80 MHz and 800 MHz the higher frequency range is applicable.

NOTE 2 These guidelines may not be applicable in all cases. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Theoretically, it is not possible to exactly predict the field strengths of fixed transmitters such as base stations for radio telephones and land mobile radios, amateur radio stations, AM and FM radio and TV broadcasting. To determine the electromagnetic environment in relation to the fixed transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location where the REHAB SERIES device is to be used exceeds the above compliance levels, the REHAB SERIES device should be monitored in order to ensure that it is functioning as intended. If unusual features are noticed, additional measures may be necessary such as re-orienting or relocating the REHAB SERIES device.

Above the frequency range from 150 kHz to 80 MHz the field strength should be less than 3 V/m.

Recommended separation distances between portable and mobile RF telecommunications 4- series and the REHAB SERIES device

The REHAB SERIES device is intended for operation in an electromagnetic environment where RF disturbances are monitored. The customer or user of the REHAB SERIES device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF telecommunications 4- series (transmitters) and the REHAB SERIES device – according to the output power of the communications device, as indicated below.

Rated output of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d= 1.2 \sqrt{P}$	80 MHz to 800 MHz $d= 0.35 \sqrt{P}$	800 MHz to 2.5 GHz $d= 0.7 \sqrt{P}$
0.01	0.12	0.035	0.07
0.1	0.38	0.11	0.22
1	1.2	0.35	0.70
10	3.8	1.1	2.2
100	12	3.5	7

For transmitters rated at a maximum output which is not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the respective column, whereby P is the maximum rated output of the transmitter in Watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz the higher frequency range is applicable.

11. Ordering information

Model numbers:

11.1 CT2200

Standard Accessories:

Serial No.	Name	Quantity
1053283	CT2200 mainframe	1
1183323	Patient Interrupt Switch	1
3130805600	Mains Power cable	1
1811361	5cm ² ultrasound applicator	1
2240000006	Ultrasound Transmission Gel	1
7100100001	Rubber electrodes(60x90mm)	2
7100100000	Rubber electrodes(70x110mm)	2
9051650011	Electrode Sponges(70x100mm)	2
9051650010	Electrode Sponges(80x120mm)	2
7100000081	Self-adhesive Electrodes(50x50mm)	4
7100000152	Self-adhesive Electrodes(50x100mm)	4
7200300010	Fixation strap(75x1200mm)	1
7200300050	Fixation strap(75x600mm)	1
7101000016	Stim Lead Wires	2
2201132900	User manual	1
1181149000	Combo Lead Wires	1

Optional Accessories:

Serial No.	Name	Quantity
1811373	1cm ² ultrasound applicator	1

11.2 MT2200

Standard Accessories:

Serial No.	Name	Quantity
1093294	MT2200 mainframe	1
1183323	Patient Interrupt Switch	1
3130805600	Mains Power cable	1
7100100001	Rubber electrodes(60x90mm)	2
7100100000	Rubber electrodes(70x110mm)	2
9051650011	Electrode Sponges(70x100mm)	2
9051650010	Electrode Sponges(80x120mm)	2
7100000081	Self-adhesive Electrodes(50x50mm)	4
7100000152	Self-adhesive Electrodes(50x100mm)	4
7200300010	Fixation strap(75x1200mm)	1
7200300050	Fixation strap(75x600mm)	1
7101000016	Stim Lead Wires	2
2201132900	User manual	1

11.3 UT2200

Standard Accessories:

Serial No.	Name	Quantity
1043305	UT2200 mainframe	1
3130805600	Mains Power cable	1
1811361	5cm ² ultrasound applicator	1
2240000006	Ultrasound Transmission Gel	1
2201132900	User manual	1

11.4 MTM200

Standard Accessories

Serial No.	Name	Quantity
1223317	MTM200 mainframe	1
1183323	Patient Interrupt Switch	1
7100100001	Rubber electrodes(60x90mm)	2
7100100000	Rubber electrodes(70x110mm)	2
9051650011	Electrode Sponges(70x100mm)	2
9051650010	Electrode Sponges(80x120mm)	2
7100000081	Self-adhesive Electrodes(50x50mm)	4
7100000152	Self-adhesive Electrodes(50x100mm)	4
7200300010	Fixation strap(75x1200mm)	1
7200300050	Fixation strap(75x600mm)	1
7101000016	Stim Lead Wires	2

11.5 BTM200

Standard Accessories:

Serial No.	Name	Quantity
1223315	BTM200 mainframe	1

11.6 VAM200

Standard Accessories:

Serial No.	Name	Quantity
1223314	VAM200 mainframe	1
7900033220	Vacuum Electrode Cups Ø 60mm	4
7120033290	Vacuum Sponges Ø 62mm	4
7130033300	Vacuum Lead Hose (Red)	2
7130033310	Vacuum Lead Hose (Black)	2

11.7 EMG200

Standard Accessories

Serial No.	Name	Quantity
1223316	EMG200 mainframe	1
710000081	Self-adhesive Electrodes(50x50mm)	8
10260041	Intravaginal Probe(26.5mm)	1
710100016	Stim Lead Wires	2
710100017	EMG Leadwire	1

Optional Accessories:

Serial No.	Name	Quantity
11811670	Intravaginal Probe(29.5mm)	1
10660042	Anal Probe	1

11.8 Cart:

Serial No.	Name	Quantity
1243319	Therapy System Cart	1

Manufacturer:

Name: Shenzhen Dongdixin Technology Co., Ltd.

Add: No.3 Building, Fanshen Xusheng Industrial Estate,
XilixiaoBaimang, Nanshan, Shenzhen 518108, China.

Tel: 0086-755-27652316

E-mail: service@nutekmedical.com

Authorized EC-representative:

Shanghai International Holding Corp. GmbH (Europe)

Eiffestraße 80, 20537 Hamburg Germany

Tel: 0049-40-2513175 Fax: 0049-40-255726